



Ergo Joystick Ordering Guide

Email, Fax, or call with the following information to:
Assistive Technology Engineering and Consultation Inc. (Atec)

Email Sales@ATEngineering.us

Phone 651 315-7776

Fax 612 435-2372



Aero



Original




Stingray

1) Ergojoystick Handle

<i>Part number</i>	<i>Model</i>	<i>Cost</i>
<input type="radio"/> AT00020100	Stingray	\$145
<input type="radio"/> AT00020101	Stingray with custom text	\$175
<input type="radio"/> AT00020102	Aero	\$140
<input type="radio"/> AT00020103	Original	\$130
<input type="radio"/> AT00020104	Original with custom text	\$160

2) Palm size

- _____
- Measure across the palm below the knuckles and above the thumb (see sketch).
 Write the measurement in tenths of an inch or mm (i.e. 3.3 inches or 84mm).



3) Hand

- Right
- Left

4) Joystick Shaft size

- 3/16 inch (4.8 mm) Typical for Invacare wheelchairs
- 1/4 inch (6.4 mm) Typical for Permobil, Pride, and Sunrise wheelchairs
- Other _____

When measuring the shaft size of your wheelchair controller, be sure to remove the rubber cap or plastic knob that comes from the manufacturer. You will be measuring the metal shaft underneath. For a quick check, 1/4 inch is about the diameter of a #2 pencil.



5) Color (Please Circle)

White, Black, Hot Pink, Coral Red, Royal Blue, Violet Purple

6) Custom Text

You can choose any text you want up to 2 lines of 15 characters.

- _____
- _____

Text not available on the Aero

Atec will email a quote to you and once the order is placed, the joystick will be built (by 3D printing) and shipped out in 1 month. Thank you! Please don't hesitate to contact us if you have any questions.

REVISIONS				
Rev	Description	ECO	DATE	BY
1.0	Release		3-10-2014	M Stone, Atec Inc.